




Product Datasheet

PYY Antibody (orb1264317)

Description	PYY Antibody	
Species/Host	Rabbit	 Western blot analysis of PYY Antibody in...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IF, IHC-P, WB	
Immunogen	This PYY antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 67-94 amino acids from the C-terminal region of human PYY.	
Target	PYY	 PYY Antibody IHC analysis in formalin fi...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Confocal immunofluorescent analysis of P...
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	11 kDa	
Uniprot ID	P10082	
NCBI	P10082	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	