



Product Datasheet

PHF13 Antibody (orb1264333)

Description	PHF13 Antibody	m.CHO
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This PHF13 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 263-291 amino acids from the C-terminal region of human PHF13.	
Target	PHF13	Western blot analysis of PHF13 Antibody ...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	PHF13 Antibody IHC analysis in formalin ...
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	34 kDa	
Uniprot ID	Q86YI8	
NCBI	Q86YI8	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	