




Product Datasheet

ADM Antibody (orb1264358)

Description	ADM Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, WB	
Immunogen	This ADM antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 69-96 amino acids from the Central region of human ADM.	
Target	ADM	Western Blot at 1:2000 dilution + MCF-7 ...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:2000For FACS starting dilution is: 1:25	Overlay histogram showing A549 cells sta...
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	20 kDa	
Uniprot ID	P35318	
NCBI	P35318	
Dilution Range	For WB starting dilution is: 1:2000For FACS starting dilution is: 1:25	Western Blot at 1:1000- 1:2000 dilution L...
Expiration Date	12 months from date of receipt.	