


Product Datasheet

NASP Antibody (orb1264403)


Description

NASP Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	This NASP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 152-180 amino acids from the N-terminal region of human NASP.
Target	NASP
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	85 kDa
Uniprot ID	P49321
NCBI	P49321
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western Blot
at 1:1000
dilution Lane
1: ...



Western blot
analysis of
NASP using
rabb...



Flow
cytometric
analysis of
HL-60 cells ...