

Product Datasheet

HIST1H2AG Antibody (orb1264408)

Description
HIST1H2AG Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This HIST1H2AL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 103-130 amino acids from the C-terminal region of human HIST1H2AL.

Target

HIST1H2AG

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


14 kDa

Uniprot ID
P0C0S8
NCBI
P0C0S8
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in
HL-60,
HepG2, K...



HIST1H2AL
Antibody IHC
analysis in
forma...



Flow
cytometric
analysis of
HepG2 cells

...