

Product Datasheet

PLA2G1B Antibody (orb1264414)

Description

PLA2G1B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This PLA2G1B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 116-145 amino acids from the C-terminal region of human PLA2G1B.

Target

PLA2G1B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

16 kDa

Uniprot ID
[P04054](#)
NCBI
[P04054](#)
Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

m.spleen


55



Western blot analysis of PLA2G1B Antibod...



Western blot analysis of PLA2G1B using r...



Flow cytometric analysis of 293 cells (r...