



Product Datasheet

TSPAN7 Antibody (orb1264416)

Description	TSPAN7 Antibody	
Species/Host	Rabbit	 Western blot analysis of TSPAN7 Antibody...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, WB	
Immunogen	This TSPAN7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 98-124 amino acids from the Central region of human TSPAN7.	
Target	TSPAN7	 Flow cytometric analysis of CEM cells (r...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	28 kDa	
Uniprot ID	P41732	
NCBI	P41732	
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	