

## Product Datasheet

### TXNDC12 Antibody (orb1264417)

## Description

TXNDC12 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This TXNDC12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 144-172 amino acids from the C-terminal region of human TXNDC12.

### Target

TXNDC12

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

19 kDa

### Uniprot ID

[O95881](#)

### NCBI


[O95881](#)

### Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

### Expiration Date


12 months from date of receipt.



Western blot analysis of TXNDC12 Antibod...



Western blot analysis of TXNDC12 Antibod...



TXNDC12 Antibody IHC analysis in formali...