



Product Datasheet

CWC15 Antibody (orb1264435)

Description	CWC15 Antibody	
Species/Host	Rabbit	
Reactivity	Hamster, Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	Western blot analysis of CWC15 Antibody ...
Immunogen	This CWC15 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 123-152 amino acids from the Central region of human CWC15.	
Target	CWC15	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	CWC15 Antibody IHC analysis in formalin ...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	27 kDa	
Uniprot ID	Q9P013	
NCBI	Q9P013	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	