

Product Datasheet

CTC1 Antibody (orb1264456)

Description

CTC1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This CQ068 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 767-795 amino acids from the Central region of human CQ068.

Target

CTC1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


135 kDa

Uniprot ID
Q2NKJ3
NCBI
Q2NKJ3
Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis of CQ068 Antibody ...



Flow cytometric analysis of K562 cells (...)