

Product Datasheet

TFRC Antibody (orb1264489)

Description

TFRC Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This CD71 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 649-677 amino acids from the C-terminal region of human CD71.

Target

TFRC

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


85 kDa

Uniprot ID
P02786
NCBI
P02786
Dilution Range


For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Antibody staining CD71 in human placenta...



Western Blot at 1:2000 dilution Lane 1: ...



Western blot analysis of CD71 Antibody i...