




Product Datasheet

IMMT Antibody (orb1264500)

Description	IMMT Antibody	
Species/Host	Rabbit	 Western blot analysis of IMMT Antibody i...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This IMMT antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 700-729 amino acids from the C-terminal region of human IMMT.	
Target	IMMT	 Western blot analysis of IMMT using rabb...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Flow cytometric analysis of hela cells u...
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	84 kDa	
Uniprot ID	Q16891	
NCBI	Q16891	
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	