

Product Datasheet

DLGAP5 Antibody (orb1264587)

Description

DLGAP5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This DLG7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 228-258 amino acids from the Central region of human DLG7.

Target

DLGAP5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

95 kDa

Uniprot ID
[Q15398](#)
NCBI
[Q15398](#)
Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


250

 130
95
72
55

Western blot analysis of DLG7 Antibody i...



DLG7 Antibody IHC analysis in formalin f...



Flow cytometric analysis of 293 cells (r...