

Product Datasheet

M6PR Antibody (orb1264599)

Description

M6PR Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, WB

Applications
Immunogen

This M6PR antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 230-260 amino acids from the C-terminal region of human M6PR.

Target

M6PR

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

31 kDa

Uniprot ID
P20645
NCBI
P20645
Dilution Range

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


130

55
36
28

Western blot analysis in A549 cell line ...



Flow cytometric analysis of A549 cells (...)



Confocal immunofluorescent analysis of M...