

## Product Datasheet

### COL9A1 Antibody (orb1264616)

## Description

COL9A1 Antibody

## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

IHC-P, WB

## Applications

## Immunogen

This COL9A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 428-456 amino acids from the Central region of human COL9A1.

## Target

COL9A1

## Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

## Form/Appearance

Liquid

## Concentration

batch dependent

## Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

## Note

For research use only

## Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

## Isotype

Rabbit Ig

## Clonality

Polyclonal

## MW

92 kDa

## Uniprot ID

[P20849](#)

## NCBI


[P20849](#)

## Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

## Expiration Date


12 months from date of receipt.



Western Blot  
at 1:1000  
dilution Lane  
1: ...



Western blot  
analysis in  
MCF-7 cell  
line...



Formalin-  
fixed and  
paraffin-  
embedded  
hum...