

Product Datasheet

LECT1 Antibody (orb1264644)

Description

LECT1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This LECT1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 253-281 amino acids from the C-terminal region of human LECT1.

Target

LECT1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



37 kDa

Uniprot ID
075829
NCBI
075829
Dilution Range


For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Western blot analysis of anti-LECT1 Anti...



Western blot analysis of LECT1 using rab...