

Product Datasheet

IDO1 Antibody (orb1264645)

Description

IDO1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This INDO antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 79-105 amino acids from the Central region of human INDO.

Target

IDO1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

45 kDa

Uniprot ID

P14902

NCBI

P14902

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

130

55

36

28

Western blot analysis
in mouse cerebellu...

95


72

55

36

28

Western blot analysis
in K562 cell line ...



INDO antibody
immunohistochemistry
analy...