


Product Datasheet

Jph3 Antibody (orb1264659)

Description	Jph3 Antibody	
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This Junctophilin 3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 549-578 amino acids from the Central region of human Junctophilin 3.	Western blot analysis of anti-Junctophilin...
Target	Jph3	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Formalin-fixed and paraffin-embedded hum...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	81 kDa	
Uniprot ID	Q9ET77	
NCBI	Q9ET77	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	