



Product Datasheet

MSX2 Antibody (orb1264669)

Description	MSX2 Antibody	
Species/Host	Rabbit	 Western blot analysis of Msx2 using rabb...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This Msx2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 191-222 amino acids from the C-terminal region of human Msx2.	
Target	MSX2	 Formalin-fixed and paraffin-embedded hum...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	29 kDa	
Uniprot ID	P35548	
NCBI	P35548	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	