

Product Datasheet

CCND1 Antibody (orb1264674)

Description

CCND1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This Cyclin D1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 68-97 amino acids from human Cyclin D1.

Target

CCND1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

[P24385](#)

NCBI


[P24385](#)

Dilution Range


For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Antibody staining
Cyclin D1 in
human ton...



Western Blot
at 1:2000
dilution Lane
1: ...



Western Blot
at 1:2000
dilution Lane
1: ...