





Product Datasheet

ACOX1 Antibody (orb1264692)

Description	ACOX1 Antibody	
Species/Host	Rabbit	 Western blot analysis of ACOX1 using rab...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This ACOX1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 629-659 amino acids from the C-terminal region of human ACOX1.	 Western blot analysis of ACOX1 using rab...
Target	ACOX1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	 Western blot analysis of ACOX1 using rab...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	74 kDa	
Uniprot ID	Q15067	
NCBI	Q15067	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	