



Product Datasheet

ZMPSTE24 Antibody (orb1264717)

Description	ZMPSTE24 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This ZMPSTE24 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 440-470 amino acids from the C-terminal region of human ZMPSTE24.	Antibody is used in Western blot to detect...
Target	ZMPSTE24	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Formalin-fixed and paraffin-embedded hum...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	55 kDa	
Uniprot ID	O75844	
NCBI	O75844	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	