

Product Datasheet

SLC2A8 Antibody (orb1264757)

Description

SLC2A8 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This SLC2A8 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 260-292 amino acids from the Central region of human SLC2A8.

Target

SLC2A8

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

51 kDa

Uniprot ID

[Q9NY64](#)

NCBI

[Q9NY64](#)


Dilution Range

For WB starting dilution is: 1:2000

Expiration Date

12 months from date of receipt.

Western Blot



Western Blot
at 1:2000
dilution Lane
1: ...