


Product Datasheet

ATP1A1 Antibody (orb1264762)

Description	ATP1A1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This NaK ATPase antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 12-46 amino acids from human NaK ATPase.
Target	ATP1A1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:2000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	113 kDa
Uniprot ID	P05023
NCBI	P05023
Dilution Range	For WB starting dilution is: 1:2000
Expiration Date	12 months from date of receipt.



Western Blot
at 1:2000
dilution +
293 wh...