

Product Datasheet

Cox4i1 Antibody (orb1264779)

Description

Cox4i1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This COXIV antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 115-169 amino acids from the mouse region of human COXIV.

Target

Cox4i1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

20 kDa

Uniprot ID

[P19783](#)

NCBI


[P19783](#)

Dilution Range

For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000


Expiration Date

12 months from date of receipt.



Overlay histogram showing NIH/3T3 cells

...



Western Blot at 1:2000 dilution Lane 1: ...