

Product Datasheet

ATP5G1 Antibody (orb1264786)

Description

ATP5G1 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This ATP5G1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 41-71 amino acids from the Central region of human ATP5G1.

Target

ATP5G1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

14 kDa

Uniprot ID

[P05496](#)

NCBI


[P05496](#)

Dilution Range

For WB starting dilution is: 1:2000

Expiration Date

12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...