


Product Datasheet

PCDH10 Antibody (orb1264907)

Description	PCDH10 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This PCDH10 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 128-160 amino acids from human PCDH10.
Target	PCDH10
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000-1:2000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	113 kDa
Uniprot ID	Q9P2E7
NCBI	Q9P2E7
Dilution Range	For WB starting dilution is: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western Blot
at 1:1000-
1:2000
dilution L...