

Product Datasheet

KCTD11 Antibody (orb1264927)

Description

KCTD11 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This KCTD11 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 21-53 amino acids from human KCTD11.

Target

KCTD11

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25 For WB starting dilution is: 1:1000-1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

26 kDa

Uniprot ID

[Q693B1](#)

NCBI


[Q693B1](#)

Dilution Range


For IHC-P starting dilution is: 1:25 For WB starting dilution is: 1:1000-1:2000

Expiration Date

12 months from date of receipt.



Antibody staining
KCTD11 in
human
kidney...



Western Blot
at 1:1000-
1:2000
dilution L...