

Product Datasheet

UBE2V1 Antibody (orb1264951)

Description

UBE2V1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This UBE2V1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 113-145 amino acids from the C-terminal region of human UBE2V1.

Target

UBE2V1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

16 kDa

Uniprot ID

[Q13404](#)

NCBI

[Q13404](#)


Dilution Range

For WB starting dilution is: 1:1000


Expiration Date

12 months from date of receipt.

m.cerebellum
72



Western blot analysis of UBE2V1 Antibody...



Western blot analysis of UBE2V1 using ra...