



Product Datasheet

USP3 Antibody (orb1264973)

Description	USP3 Antibody	
Species/Host	Rabbit	 Western blot analysis of hUSP3-Y505 in K...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This USP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 489-519 amino acids from the C-terminal region of human USP3.	
Target	USP3	 Formalin-fixed and paraffin-embedded hum...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	59 kDa	
Uniprot ID	Q9Y6I4	
NCBI	Q9Y6I4	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	