


Product Datasheet

GDF11 Antibody (orb1265001)


Description

GDF11 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This GDF11 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 32-61 amino acids from the N-terminal region of human GDF11.
Target	GDF11
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	O95390
NCBI	O95390
Dilution Range	For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Formalin-
fixed and
paraffin-
embedded
hum...



Formalin-
fixed and
paraffin-
embedded
hum...