

Product Datasheet

LYZ Antibody (orb1265006)

Description

LYZ Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This LYZ antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 119-154 amino acids from the C-terminal region of human LYZ.

Target

LYZ

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:25

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

17 kDa

Uniprot ID

[P61626](#)

NCBI


[P61626](#)

Dilution Range


For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:25

Expiration Date


12 months from date of receipt.



Western Blot at
1:2000 dilution Lane
1: ...



Immunohistochemical
analysis of paraffin...



Immunohistochemical
analysis of paraffin...