

Product Datasheet

GADD45GIP1 Antibody (orb1265024)

Description

GADD45GIP1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This GADD45GIP1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 13-45 amino acids from the N-terminal region of human GADD45GIP1.

Target

GADD45GIP1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

25 kDa

Uniprot ID

[Q8TAE8](#)

NCBI


[Q8TAE8](#)

Dilution Range

For FACS starting dilution is: 1:25For WB starting dilution is: 1:1000


Expiration Date

12 months from date of receipt.



Flow
cytometric
analysis of
MCF-7 cells

...



Western blot
analysis of
lysate from
293...