

Product Datasheet

MPP7 Antibody (orb1265027)

Description

MPP7 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This MPP7 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 2-36 amino acids from the N-terminal region of human MPP7.

Target

MPP7

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For FACS starting dilution is: 1:25For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

66 kDa

Uniprot ID

[Q5T2T1](#)

NCBI


[Q5T2T1](#)

Dilution Range


For IHC-P starting dilution is: 1:25For FACS starting dilution is: 1:25For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Flow cytometric analysis of Hela cells
u...



Immunohistochemical analysis of paraffin...