

Product Datasheet

ECE2 Antibody (orb1265070)

Description

ECE2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This ECE2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 551-579 amino acids from the Central region of human ECE2.

Target

ECE2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

100 kDa

Uniprot ID

[O60344](#)

NCBI


[O60344](#)

Dilution Range


For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis in U251 cell line ...