

Product Datasheet

MAP2 Antibody (orb1265145)

Description

MAP2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This MAP2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1766-1796 amino acids from the C-terminal region of human MAP2.

Target

MAP2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

200 kDa


Uniprot ID
P11137
NCBI
P11137
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100


Expiration Date

12 months from date of receipt.

250 1 2



Western blot analysis of hMAP2-Q425 using...



Formalin-fixed and paraffin-embedded hum...