


Product Datasheet

FUNDC2 Antibody (orb1265161)

Description	FUNDC2 Antibody	 Western blot analysis in mouse kidney ti...
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This FUNDC2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 122-150 amino acids from the C-terminal region of human FUNDC2.	
Target	FUNDC2	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	21 kDa	
Uniprot ID	Q9BWH2	
NCBI	Q9BWH2	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	