


Product Datasheet

Txnip Antibody (orb1265174)

Description	Txnip Antibody	m.stomach
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This Mouse Txnip antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 116-145 amino acids from the Central region of mouse Txnip.	
Target	Txnip	Western blot analysis in mouse stomach t...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	44 kDa	
Uniprot ID	Q8BG60	
NCBI	Q8BG60	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	