

Product Datasheet

POLDIP3 Antibody (orb1265183)

Description

POLDIP3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This POLDIP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 93-120 amino acids from the N-terminal region of human POLDIP3.

Target

POLDIP3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For IF starting dilution is: 1:25For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

46 kDa

Uniprot ID

[Q9BY77](#)

NCBI


[Q9BY77](#)

Dilution Range


For IHC-P starting dilution is: 1:25For IF starting dilution is: 1:25For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Fluorescent image of HeLa cells stained ...



Western blot analysis in Jurkat, mouse C...