

Product Datasheet

ATP1A3 Antibody (orb1265232)

Description

ATP1A3 Antibody

NCI-H292

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This ATP1A3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 805-833 amino acids from the Central region of human ATP1A3.

Target

ATP1A3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

112 kDa

Uniprot ID

[P13637](#)

NCBI


[P13637](#)

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in NCI-H292 cell l...