

Product Datasheet

LEO1 Antibody (orb1265345)

Description

LEO1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This LEO1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 130-159 amino acids from the N-terminal region of human LEO1.

Target

LEO1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000-1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

75 kDa

Uniprot ID

[Q8WVC0](#)

NCBI

[Q8WVC0](#)

Dilution Range


For WB starting dilution is: 1:1000-1:2000

Expiration Date


12 months from date of receipt.



Western Blot
at 1:1000-
1:2000
dilution L...



Western blot
analysis of
anti-Leo1
Antib...



Western blot
analysis of
Leo1 using
rabb...