


Product Datasheet

A20A1 antibody (orb126535)

Description	Rabbit Polyclonal to A20A1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q5TYW2
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 3...