

## Product Datasheet

### LGALS3BP Antibody (orb1265350)

## Description

LGALS3BP Antibody

## Species/Host

Rabbit

## Reactivity

Mouse

## Conjugation

Unconjugated

## Tested

WB

## Applications

## Immunogen

This LGALS3BP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 383-412 amino acids from the Central region of human LGALS3BP.

## Target

LGALS3BP

## Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

## Form/Appearance

Liquid

## Concentration

batch dependent

## Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

## Note

For research use only

## Application notes

For WB starting dilution is: 1:1000

## Isotype

Rabbit Ig

## Clonality

Polyclonal

## MW

65 kDa

## Uniprot ID

[Q08380](#)

## NCBI


[Q08380](#)

## Dilution Range

For WB starting dilution is: 1:1000

## Expiration Date

12 months from date of receipt.



Western blot analysis in mouse lung tiss...