


Product Datasheet

ZNF185 Antibody (orb1265407)

Description	ZNF185 Antibody	 Western blot analysis in NCI-H460 cell l...	
Species/Host	Rabbit		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	WB		
Immunogen	This ZNF185 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 155-183 amino acids from the N-terminal region of human ZNF185.		
Target	ZNF185		
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.		
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		
Note	For research use only		
Application notes	For WB starting dilution is: 1:1000		
Isotype	Rabbit Ig		
Clonality	Polyclonal		
MW	74 kDa		
Uniprot ID	O15231		
NCBI	O15231		
Dilution Range	For WB starting dilution is: 1:1000		
Expiration Date	12 months from date of receipt.		