

## Product Datasheet

### ST6GALNAC5 Antibody (orb1265500)

## Description

ST6GALNAC5 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

This ST6GALNAC5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 110-138 amino acids from the C-terminal region of human ST6GALNAC5.

### Target

ST6GALNAC5

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

38 kDa

### Uniprot ID

[Q9BVH7](#)

### NCBI


[Q9BVH7](#)

### Dilution Range

For WB starting dilution is: 1:1000

### Expiration Date

12 months from date of receipt.



Western blot analysis in ZR-75-1 cell li...