


Product Datasheet

Collagen III (Cleaved-Gly1221) antibody (orb126562)

Description	Rabbit polyclonal antibody to Collagen III (Cleaved-Gly1221)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P02461
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from A...