

## Product Datasheet

### RUVBL1 Antibody (orb1265658)

## Description

RUVBL1 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

FC, IHC-P, WB

### Applications

### Immunogen

This RUVBL1 (TIP49) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 247-280 amino acids from the Central region of human RUVBL1 (TIP49).

### Target

RUVBL1

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

50 kDa

### Uniprot ID

[Q9Y265](#)

### NCBI


[Q9Y265](#)

### Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

### Expiration Date


12 months from date of receipt.



Western blot analysis of RUVBL1 using ra...



RUVBL1 (TIP49) Antibody immunohistochemi...



Flow cytometric analysis of 293 cells (r...