

Product Datasheet

C1R (L chain, Cleaved-Ile464) antibody (orb126570)

Description Rabbit Polyclonal to C1R (L chain, Cleaved-Ile464).

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Internal

Concentration 1mg/ml


Note For research use only

Clonality Polyclonal

Uniprot ID **P00736**

Dilution Range WB: 1:500-1:3000, ELISA: 1:20000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts from 2...