


Product Datasheet

Collagen IV (Cleaved-Leu1425) antibody (orb126576)

Description	Rabbit polyclonal antibody to Collagen IV
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q01955
Dilution Range	WB: 1:500-1:3000, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western
blot
analysis
of
extracts
from 293T cells