

Product Datasheet

RNASE2 Antibody (orb1265800)

Description

RNASE2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This RNASE2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 124-152 amino acids from the C-terminal region of human RNASE2.

Target

RNASE2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

18 kDa

Uniprot ID

P10153

NCBI


P10153

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of lysate from hum...